

BP37™

Fitting Instructions BABYPOD - Variant



IMPORTANT:
TO USE THE BP37 WITH BABYPODS,
REFER TO APPROPRIATE IFU'S.




figure 1

Locate Controller
mounting clip and
connect with metal clip.




figure 2




figure 3




figure 4

Rotate Controller as indicated when
attaching Controller. (figures 2 & 3)

Ensure cable is not trapped
between Controller and
BabyPod shell.



Check Controller is connected securely.




figure 5

To remove the Controller, you lift the
Controller in a circular motion and slide
the Controller away from the BabyPod.



Manufactured in France on behalf of:



Advanced Healthcare Technology Ltd.
Theobalds Enterprise Centre, Platinum Way, Waltham Cross, Herts., EN8 8YD, United Kingdom
Tel: +44(0)1992 535933 Email: sales@babypod.com
www.babypod.com



Advanced Healthcare Technologies s.r.o.
Varšavská 715/36, 120 00 Prague 2, Czech Republic



BP37™ is a patented product.
BP37™ is a trademark of Advanced Healthcare Technology Ltd.