

BP37™

Fitting Instructions

LS - Variant



IMPORTANT:
TO USE THE BP37,
REFER TO APPROPRIATE IFU'S.




Parts Supplied:

Rail Clamp (Part No. BP37-014)
2 x M5 x12mm screws


Tools Required (not supplied):

4mm Allen Key




Attach Rail Clamp to Controller
with supplied Screws.

Turn Rail Clamp handle 90° to secure
Controller to Medi-Rail.



Use of Cable guide



Push cable guide onto cable
to form a 90° angle.



Manufactured in France
on behalf of:



Advanced Healthcare Technology Ltd.
Theobalds Enterprise Centre,
Platinum Way, Waltham Cross,
Herts., EN8 8YD, United Kingdom
Tel: +44(0)1992 535933
Email: sales@babypod.com
www.babypod.com



Advanced Healthcare Technologies s.r.o.
Varšavská 715/36, 120 00 Prague 2, Czech Republic



BP37™ is a patented product.
BP37™ is a trademark of Advanced Healthcare Technology Ltd.